


▶ Basic Application Tips



End to End Application: Tear backing 2-3" from end of tape and remove from end of tape only. Apply to skin with no stretch and rub to activate adhesive. Apply center of tape with desired stretch. Finish with no stretch in final 2".



Middle Stretch Application: Tear backing across middle and begin to peel back from center. Stretch tape as indicated and apply from center towards ends. Apply final 2" of each end with no stretch.

Tennis Elbow



1.



2.

1. Extend arm, point hand down and rotate outward. Anchor tape above elbow and apply, with no stretch, to just above wrist.
2. Apply 2nd piece directly over painful area. Stretch tape 50% in middle, no stretch in ends.